



active
rehab

Rifton HTS Shower/ commode



Hygienic, simple,
affordable and
adaptable.

One seat, many possibilities

The HTS can be configured in myriad ways to fit your clients. Options range from the basic model with seat, back, base and armrests to the fully accessorized version.



Any client, any toilet, any situation

For clients, each element of the HTS creates the best possible conditions for hygiene. Its wide range of adjustment and varied support options help clients maintain positioning that promotes effective toileting, and the tilt-in-space feature and adjustable headrest make showering easier.

For caregivers, the HTS is the paragon of adaptability. It can be used on or over a round front or elongated toilet, or as a free-standing commode. It can adjust to children as they grow, and includes components and supports for clients of every size and need.



Better mobility

The large caster option improves mobility and has better resistance to picking up debris.



Accessible for clients

Features like a gas-assisted spring tilt adjustment, removable armrests, and a flip-up footboard make sit-to-stand transitions easy.



Foam molded cushions

Foam pads on the back and seat provide comfort for your clients and easy cleaning for caregivers.

The perfect shower chair

with a single product for both toileting and showering. The aluminum, stainless steel and plastic construction withstands bleach and water. The HTS tilt-in-space base tilts 5° forward or 25° back for showering and hair washing.



Grows with you

The HTS comes in three sizes, each with a wide range of adjustment. As your child grows the chair can adapt, giving you years of useful service.

Accessories



Size	Seat width	Backrest height	Distance between laterals	Tilt in space	Max user weight
Small	8-12"	16.5"	6-11.5"	-5° to 25°	34kg
Medium	10-14"	19.5"	8-13.5"		68kg
Large	13-17"	24.5"	11.5-17"		113kg