



active  
rehab

# Wonder Seat

Paediatric Seating for the Bingo OT



REHAB  
EQUIPMENT,  
TAILORED  
FOR YOUR LIFE.

# A breakthrough in rehabilitative stroller seating

A special edition of SPEX seating has been carefully developed so therapists can now provide more comfort and superior positioning for young children with complex seating needs with the freedom of a stroller rather than a wheelchair.



## SPEX POSITIONING CUSHION BASE

Strategic positioning pads can be added or removed from underneath the cushion to allow for different contours and shapes. This way we can meet the client's posture requirements present today; and readjust them in the future.



## SPEX CONTOURING SYSTEM

Pockets behind the back-foam inner create a modular and adjustable solution for changing the contour to suit different user shapes and posture needs.



## AXIAL PELVIC SUPPORTS

Multi-adjustable pelvic supports can orientate in width and angle for optimum pelvic control. These pads are mounted to the quick-release hardware meaning transfers are easy.



## AXIAL LATERAL TRUNK SUPPORTS

Comfortable lateral supports that are width, depth and multiple angle adjustable. Patented mechanically rotational system allows the lateral pad to rotate to correctly align to the user's trunk. Swing-away hardware make it easy for carers to transfer the child in and out.



## SPEX SUPERHIGH CUSHION

Specially designed for a young child who has moderate to high positioning requirements. SPEX Superhigh offers optimal adduction and abduction contours to effectively position the pelvis and legs. Spex positioning technology allows further adjustment for pelvic obliquity, hip flexion and other clinical requirements.



## FEET POSITIONING SANDALS

Straps can be fixed in a variety of ways to achieve the most appropriate support for each child. Independently adjustable to exactly accommodate the optimum foot position.



# Adjustable seating adjustable stroller

Cobra height-adjustable indoor base. The seat easily transfers into the Cobra for feeding and activity times. Tilt setting can still be adjusted.



SEAT IN FORWARD FACING DIRECTION



SEAT IN REAR FACING DIRECTION



SUN CANOPY

A smart-looking stroller that enables parents to take their child outside and explore!

## Superior comfort and style

With a wide range of coloured seating options to choose from, the SPEX Wondeseat offers a colour to suit any style. Stand out from the crowd in with our Chilli red or look pretty in pink with the Calypso.



# LIVE LIFE

PH: 1800 987 663

[support@activerehab.net.au](mailto:support@activerehab.net.au)

[www.activerehab.net.au](http://www.activerehab.net.au)

Like us on Facebook ( Life beyond the norm) 