



active
rehab

Zeal Seat

Script form 2022



**REHAB
EQUIPMENT,
TAILORED
FOR YOUR LIFE.**

Zeat Seat

Script form 2022

Client Details

Client Name

Date:

Notes

Product Details

Options Available	<input type="checkbox"/> Zeat Seat package black ZSB001	<input type="checkbox"/>
	<input type="checkbox"/> Zeat Seat package timber ZST001	<input type="checkbox"/>

Seat Width	40cm
Seat Depth	9-36cm
Seat height from floor	40-68cm
Seat angle	0, 5, 20 degrees
Backrest height	40-50cm
Backrest angle	65-115 degrees
Footplate height range	8-56cm
Footplate angle	0, 5, 20 degrees
Tray size	50x50cm
Tray height from floor	77-88cm
Overall height	up to 112cm
Overall width	55cm
Overall depth	88cm
Weight limit	75kg
width between laterals	9-24cm
Lateral size	11x12.5x1cm



Lateral Set up



Select where you want the laterals fitted

Selected where you want the laterals fitted by ticking on the slots below.

Distance you need between lateral supports

Notes for lateral supports

Zeal Seat


Script form 2022

Seat Set Up.

Select where you want the seat set up

Select how you want the seat set up by ticking in the circles and filling out the details below.

Note: the backrest mounting holes are letters, and seat and footrest mounting holes are number



Letters for backrest mounting

Numbers for seat base mounting

Numbers for footrest mounting

Numbers for tray table mounting

Notes (further information)

Zeal Seat
Brought to you by:

your logo

Active Rehab
78 Beischer St East Bendigo
Vic 3550
PH: 03 5441 8966
www.activerehab.net.au