
## marcy.

## There for you from the start.



## I am your anterior walking aid, and will assist you in your very first attempts

to walk. I can hardly wait to start off with you! I will support your balance and simultaneously train your torso muscles. You can adapt my wheelbase to suit. You can adjust my wheelbase shorter to make me more manoeuverable or by making it longer I can give you more security and support. If you find it hard to walk in a straight line, I am happy to assist you with my direction locks.

















Frame colours Size 2



Frame colours Size 3





Data	Size 0	Size 1	Size 2	Size 3
Handle height	30 - 47 cm	38 - 59 cm	53 - 75 cm	63 - 88 cm
Height of forearm rests	41 - 58 cm	49 - 70 cm	63 - 85 cm	73 - 98 cm
Handle spacing	18 - 32 cm	22 - 36 cm	30 - 44 cm	36 - 50 cm
Seat height	25 cm	30 cm	45 cm	45 cm
Total width	44 cm	48 cm	62 cm	71 cm
Total length	36.5 cm	42.5 cm	57.5 cm	61 cm
Total length (telescopic wheelbase)	39 - 47 cm	45 - 55 cm	60 - 74 cm	65 - 79 cm
Max. load	20 kg	30 kg	50 kg	65 kg
Weight	4.2 kg	5 kg	7.4 kg	9.1 kg



