



active
rehab

Klip Walker



**Keeping
you on the
move**

The klip is a lightweight posterior gait trainer for kids, teens and young adults with special need who require support and assistance with their day-to-day walking.

Open Frame

Placed behind the user, the klip is open in the walking position and allows freedom of movement and promotes walking in an upright position. The klip can also be used in the anterior position (other than the XS size). A variety of positioning accessories are available.



Accessories

The Klip's choice of accessories accommodates children with a range of supportive positioning requirements. You can choose one or combine multiple accessories for a fully supportive upright gait trainer. Installation is simple and accessories are easily added or removed as walking skills evolve.

- ✓ Available with a fold down seat
- ✓ Available with outdoor and indoor wheels
- ✓ Comes in 2 colour options (blue or red)
- ✓ Folds up

Specifications

Size	X-small	Small	Medium	Large
Width inside handles	13.5"	14"	16"	17"
Handle height	16-19.5"	31-27"	24-32"	30-37"
Outside dimensions	23x23"	24x31"	23x35"	26x37"
Max user weight	34kg	39kg	68kg	40.8kg



Pelvic Stabilizer

Provides lateral hip support, improves alignment and balance, and promotes an upright position.



Forearm Platform

Fully adjustable for comfortable weight-bearing and support from elbow to wrist.



Cushioned Sling Seat

Adjustable sling seat provides weight bearing assistances.



Tension Control

Adjusts resistance on rear wheels to control speed.



Bumper Guards

Protects walls from marks and dents and cushions the impact felt by user.