

# Baffin Trio



active  
rehab



The Baffin TRIO multifunctional device is designed for sitting and standing



### Your multi positioning device

The Baffin Trio is here to revolutionise sitting and standing by combining the two functions into the one product.

The Baffin Trio consists of a seat, back support with the 2nd spine, side supports, footrests, head support, knee pads, hip belts, a vest, an activity table, and the moveable base.

### Your multi positioning device

The Baffin features a built in spine system which can move and twist to mimic an actual spine. This means that the baffin can be adjusted and custom moulded to fit the users body shape perfectly.

Thanks to the innovative built-in spine, side supports and pelvis correction, Baffin is a perfect device for positioning. Thanks to its simplicity, it can be used at rehabilitation centres and at home to continue the therapy. Spinal support along the whole back line helps to reduce pressure on a single vertebra. Baffin helps establishes a correct sitting position.

### Key features

- Very large adjustment range - one device can accommodate a child with a height of 135 cm and an adult up to 170 cm
- Great possibilities of passive posture correction in the sitting, standing and lying position
- The possibility of setting intermediate positions in order to obtain the greatest comfort for the patient
- the possibility of standing from any position - sitting, lying and intermediate positions.
- A system of removable pads, hip supports and armrests that facilitate the transfer of the patient to the device
- Intuitive and friendly remote control

Size	Seat height	Back length	Backrest angle	Base length	Base width	Seat depth	Seat width	Weight	Max user weight
1	61cm	61cm	80°-115°	116cm	69cm	36-50cm	27-40cm	57kg	90kg
2	61cm	61/77cm	80°-115°	122cm	77cm	50-64cm	34-45cm	67kg	90kg